INSTRUCTIONS

opal[®]PRIMO[™]

Description / Indications:

Opal® Primo™ is an orthodontic light cure bonding enhancer (resin primer) that is used to bond orthodontic brackets and bondable buccal tubes to etched enamel in conjunction with light cured adhesive. Opal Primo acts as a wetting and drying agent, contains ethyl alcohol solvent (drying agent), and contains 7.5% radiopaque filler to retard evaporation and improve coverage control.

Direct Bonding Procedure:

1. Clean enamel surfaces with pumice or a similar water-based, mild abrasive.
2. Isolate the bonding field with dri-angles, cotton rolls, and cheek retractors.
3. Etch with a 35% phosphoric acid, rinse, and dry all tooth surfaces thoroughly prior to bonding with oil and moisture free air.
4. Remove lock cap from Opal Primo syringe and thread on the Opal Primo applicator tip.
5. Carefully express a small amount of the Primo bond enhancer by gently pushing the plunger while holding the syringe barrel with opposite hand.
6. If too much enhancer is expressed on the applicator tip, a clean 2 x 2 wipe or cotton roll can be used to remove the excess enhancer. One small application should easily cover an entire quadrant.
7. Apply a thin layer of Opal Primo to the surface of etched teeth in a single quadrant moving quickly from tooth to tooth.
8. Thin the enhancer over the entire buccal surface with a very short air blast from a moisture free air syringe.
9. Apply a thin layer of light cure orthodontic adhesive to bonding surface of bracket.
10. Place brackets coated with adhesive on etched and primed tooth.
11. Remove any excess bonding material from circumference of bracket base.
12. Light cure bonding material following manufacturers' recommended curing time for adhesive and bonding light.

Precautions:

1. For Professional use only.
2. Review instructions, precautions, and MSDS before beginning treatment. Use only as directed.
3. Keep products out of heat/sunlight.
4. Use quality curing lights and recommended curing times.
5. Do not use on patients with contraindications to orthodontic treatment.
6. Cure bonding resins completely, as uncured resin materials will be deficient in physical properties, stress, and wear resistance.
7. Any resin can be sensitizing to skin or can cause allergic reactions if individual is already sensitized. If seen, wash area of tissue exposed to resin completely. If dermatitis, rash, or any other allergic reaction is observed, consult a physician.
8. Replace luer lock caps when syringes are not in use.
9. When dispensing intraorally, never force syringe plunger. If there is excess resistance to flow, replace tip with a new one and verify flow before proceeding. Use only recommended tips with each material.
10. To prevent cross-contamination, tips should not be re-used. Clean and disinfect syringes between patients.

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Our syringes have an expiration date stamped on the side of the syringe consisting of one letter and three numbers.
The letter is a lot number used for manufacturing purposes and the three numbers are the expiration date.
The first two numbers are the month, and the third number is the last number of the year.

X026 = **Lot**
X026 = **MM**
X026 = **YYYY**

EN - Refrigerate
DE - Kühlen
FR - Réfrigérer
NL - In de koelkast bewaren
IT - Conservare in frigorifero
ES - Refrigerar
PT - Refrigerar
SV - Kyla ned
DA - Opbevares i køleskab
FI - Jäähdytetävä
EL - Το διατηρείτε στο ψυγείο

EN - Do not reuse to avoid cross-contamination.
DE - Zur Vermeidung von Kreuzkontamination nicht wiederverwenden.
FR - Ne pas réutiliser pour éviter la contamination croisée.
NL - Niet opnieuw gebruiken om kruisbesmetting te vermijden.
IT - Non riutilizzare per evitare la contaminazione crociata.
ES - No reutilice para evitar la contaminación cruzada.
PT - Não reutilizar para evitar a contaminação cruzada.
SV - För att undvika korskontamination är återanvändning inte tillåten.
DA - Må ikke genbruges for at undgå krydskontaminering.
FI - Älä käytä uudelleen ristikontaminaation riskin vuoksi.
EL - Μην το επαναχρησιμοποιείτε προς αποφυγή διασταυρωτικής μόλυνσης.

EN - For prescription use only
DE - Nur auf Rezept erhältlich
FR - Pour utilisation sur ordonnance uniquement.
NL - Uitsluitend voor gebruik op recept
IT - Prescritto dal dentista.
ES - Para uso exclusivo con receta.
PT - Destina-se apenas à utilização mediante prescrição
SV - Får endast utlämmas mot recept
DA - Administreret af tandlæge
FI - Reseptivalmistie
EL - Για χρήση μόνο με συνταγή γιατρού

EN - Manufactured by
DE - Hersteller
FR - Fabriqué par
NL - Geproduceerd door
IT - Prodotto da
ES - Fabricado por
PT - Fabricado por
SV - Tillverkad av
DA - Fremstillet af
FI - Valmistaja
EL - Κατασκευάστηκε από την

EN - See Instructions
DE - Siehe Anleitung
FR - Voir le mode d'emploi
NL - Zie instructies
IT - Consultare le istruzioni
ES - Ver instrucciones
PT - Ver instruções
SV - Se instruktioner
DA - Se brugsanvisningen
FI - Lue ohjeet
EL - Δείτε τις οδηγίες

Xi

EN - Irritant
DE - Reizend
FR - Irritant
NL - Irriterend
IT - Irritante
ES - Irritante
PT - Irritante
SV - Irriterande
DA - Lokalirriterende
FI - Ärsyttävä
EL - Ερεθιστικό

F

EN - Highly Flammable
DE - Leichtentzündlich
FR - Facilement inflammable
NL - Licht ontvlambaar
IT - Facilmente infiammabile
ES - Fácilmente inflamable
PT - Facilmente inflamável
SV - Högt Brandfarligt
DA - Letantændelig
FI - Helposti Syttävä
EL - Εύφλεκτο

Manufactured by:

Opal Orthodontics

505 West 10200 South
South Jordan, Utah 84095 USA
Toll Free: 1 (888) 863-5883
Outside U.S. Call: +1(801) 553-4177

Made in USA.
www.opalorthodontics.com
U.S. Pat. No.: 5,708,052
© 2010 Opal Orthodontics. All rights reserved.

28058.3

010510